

Headspace Vials

For Headspace Autosampler Instruments



Precision screw thread ND18 Vials

Headspace Autosampler Vials

Precision Screw thread ND18 Vials

- Screw Headspace caps, septa, and vials are designed especially for use with Agilent and CTC rotating or robotic arm samplers.
- Certified by GC/MS to be contaminate free for VOC analysis.
- 10ml, 20ml Clear or Amber precision thread vial.
- Magnetic screw caps 8mm center hole.
- Septa are blue PTFE sample side; White Silicone outside.
- Septa diameter : 17.5 x 1.5mm
- Vial 10ml diameter : 22.5 x 46 mm, Vial 20ml diameter : 22.5 x 75.5mm
- Expansion 33, Clear Borosilicate Type I
- Schott Glass, Amber Borosilicate Type I
- Precision Screw thread Headspace Vials rounded bottom

Ordering Information :

Precision Screw thread ND18 Vials

Description	Part No.	Qty
10ml Screw thread ND18 Vial, Clear Round Bottom	121801	100/pk
20ml Screw thread ND18 Vial, Clear Round Bottom	121802	100/pk
20ml Screw thread ND18 Vial, Amber Round Bottom	121806	100/pk

Magnetic Screw Cap and Septa for Precision Thread Vials ND18

Description	Part No.	Qty
18mm Magnetic screw cap, silver, 8mm center hole with blue PTFE/white Silicone septa	121803	100/pk
Septa Blue PTFE/white Silicone, 17.5 x 1.5mm	121804	100/pk
18mm Magnetic screw cap, silver, 8mm centre hole	121805	100/pk

