

Nitrogen Generator Model LGN 32 MA

Flow rate : 32 liter/min



Membrane Nitrogen Generator

The LABGEN Nitrogen Generator Model LGN 32 MA is a self-contained membrane nitrogen generator that produces pure grade nitrogen at output pressure 100 psig and flow rate up to 32 liter/min. Nitrogen is produced by utilizing a combination of inbuilt silent oil free air compressor, coalescing filter and nitrogen membrane separation technology components.

Applications include: LC/MS/MS Mass Spectrometer, Nebulizer gases, chemical and solvent evaporation, instrument supply and purge, evaporative light scattering equipment.

Features & Benefits:

- No external compressor required.
- Inbuilt Silent oil free air compressor, Includes 2 years compressor warranty.
- LED Indicator light when air compressors are in operated, run time over and time reach to maintenance.
- Quiet, Insulated compressor compartment for minimum disruption.
- Nitrogen Flow capacity up to 32 Liter/min.
- Flexibility to position generator where required, the generator is supplied with caster wheels.
- Designed as plug & play system.
- Cost effective than any other gas supply method.
- Gas on demand, no health hazards, no need to worry about running out of gas.

Ordering Information :

Part No.	Descriptions
321601	Nitrogen Generator Model LGN 32 MA, 220V, 50 Hz

Spare Parts

Part No.	Descriptions
321605	Silent Oil free Air compressor assembly
321606	Intake Filter Element
321607	Element Filter for AMG Water Separator filter
321608	Element Filter for AMD Micro Mist Separator filter, 0.01 um
321609	Element Filter for AM Micro Mist Separator filter, 0.3 um

Technical Specifications

Technology	Nitrogen Membrane Separation
Inbuilt Air Compressor	Silent Oil free Air Compressor,
Nitrogen membrane Outlet Purity	95% - 99.9% (depend on flow)
Particles	< 0.01µm
Phthalates	None
Suspended Liquids	None
Noise Level	≤ 54 dB(A) @ 1m
Nitrogen Storage tank	6 gallons (22 Liters)
Gas Outlet Connection	1x 1/4" NPT Female
Electrical Requirements	220 V 50 Hz, 8.0 A
Dimensions (mm) W x D x H	600 x 760 x 720 mm
Weight (kg/ lbs)	95 / 209

Standard Accessories : 2 AMG Water Separator filter, 1 AM Micro Mist Separator Air filter 0.3um, 1 AMD Micro Mist Separator Air filter 0.01um

Spare Parts and Service Interval Component



P/N : 321605
Oil free Air compressor
820 W, 220 V/ 50Hz
Warranty 2 years



P/N : 321606
Intake Element Air Filter
12 months



P/N :321607
Element filter for AMG Water Separator
12 months



P/N : 321608
Element filter for AMD Micro Mist Separator, 0.01um
12 months



P/N : 321609
Element filter for AM Micro Mist Separator, 0.3um
12 months