

Multimodule Vortex Mixer Model LC-Vortex MT



Multimodule Vortex Mixer

Features

- The large LCD screen displays the speed and time, which is convenient for setting and reading;
- Suitable for short time (jog) or long time continuous work;
- Brushless DC motor, long life, strong power, maintenance-free;
- Special vibrating head, strong toughness, sufficient elasticity, and long service life;
- Wide speed range, the maximum speed is 3000 rpm, stepless speed regulation;
- The module is easy to install, and the eccentric bearing design is durable;
- Using high-quality foot pads, the instrument has good stability and low operating noise.

Optional (Order Separate)



Module 1
8 x 17 mm
For 20 – 25ml test tube



Module 2
8 x 13 mm
For 7 – 15ml test tube



Module 3
16 x 9 mm
For 1.5 – 5ml test tube



Module 4
15 x 7.5 mm
For 0.5ml test tube



Module 5
48 x 4 mm
For 0.2ml test tube

Specifications

Model	LC-Vortex MT
Voltage	100 - 240 VAC / 12VDC
Power	20 W
Display	LCD Screen
Oscillation Mode	Circle
Turnover Of Diameter	4.5mm
Speed Range	100-3000rpm
Display Of Speed	LCD Screen
Timing Range	0-99min
Mode Of Operation	Continuous Operation/Point Operation
Net Weight kg	2.6
Dimension mm	145 × 160 × 115

