



2019-2020

VREZ

Filtration & Separation Solutions

Filtration Manifold

Focus on
laboratory sample filtration solutions



บริษัท ซี.อี. อินสตรูमेंท์ (ประเทศไทย) จำกัด
68/876 หมู่ 8 ถนนรัตนาธิเบศร์ ต. บางกระสอ
อ. เมืองนนทบุรี จ. นนทบุรี 11000

โทรศัพท์ 02 9658971-3 โทรสาร 02 9658974


Email : sales@cethailand.com www.cethailand.com

TIANJIN SCICHRO TECHNOLOGIES CO., LTD.



FILTRATION MANIFOLD

Enjoy the Laboratory



GUIDANCE

1. choose a suitable manifold rack



2. choose a suitable filtration holder



Stainless steel
300ml filter holder
clamp-free



Stainless steel
150ml filter holder
clamp-free



Boro.glass
300ml filter holder
with clamp



PP funnel, sterile
250ml filter holder



PP funnel, sterile
250ml filter holder
with membrane filter



FEATURE

Main parts is made of SUS316

Each parts can be simple and rapidly assembled, which offer excellent chemical resistance and can be fastly sterilized by flame, steam, oven etc..

Individual control valve

Each funnel valve can be independently controlled.

Compatible to different size of holder


Manifold Branch Ports are available to fit 150ml, 300ml, 250ml, 300ml holders.

Meet various laboratory need



PRODUCT INFORMATION


Filtration manifold with 300ml SS holder



SPECIFICATION

Item	VREZ-1CF	VREZ-3CF	VREZ-6CF
Part.No.	267102-71	267303-72	267603-73
Branch number	1 Port	3 Ports	6 Ports
Holder number	1 Set	3 Sets	6 Sets
Funnel volume	300ml	300ml	300ml
Filter diameter	47-50mm	47-50mm	47-50mm
outlet	ø 8 mm	ø 8 mm	ø 8 mm
Parts	Material	Material	Material
Rack	SUS316	SUS316	SUS316
Handle	Anodized Alum.	Anodized Alum.	Anodized Alum.
Funnel & Lid	SUS316	SUS316	SUS316
Funnel support base	SUS316	SUS316	SUS316

Filtration manifold with 150ml SS holder




SPECIFICATION

Item	VREZ-1DF	VREZ-3DF	VREZ-6DF
Part.No.	267104-81	267304-82	267604-83
Branch number	1 Port	3 Ports	6 Ports
Holder number	1 Set	3 Sets	6 Sets
Funnel volume	150ml	150ml	150ml
Filter diameter	47-50mm	47-50mm	47-50mm
outlet	ø 8 mm	ø 8 mm	ø 8 mm
Parts	Material	Material	Material
Rack	SUS316	SUS316	SUS316
Handle	Anodized Alum.	Anodized Alum.	Anodized Alum.
Funnel & Lid	SUS316	SUS316	SUS316
Funnel support base	SUS316	SUS316	SUS316




Filtration manifold with 300ml Boro.glass holder



SPECIFICATION

Item	VREZ-1AF	VREZ-3AF	VREZ-6AF
Part.No.	267100-51	267300-52	267600-53
Branch number	1 Port	3 Ports	6 Ports
Holder number	1 Set	3 Sets	6 Sets
Funnel volume	300ml	300ml	300ml
Filter diameter	47-50mm	47-50mm	47-50mm
outlet	ø 8 mm	ø 8 mm	ø 8 mm
Parts	Material	Material	Material
Rack	SUS316	SUS316	SUS316
Handle	SUS316	SUS316	SUS316
Funnel & Lid	Boro.glass / SS	Boro.glass / SS	Boro.glass / SS
Funnel support base	Boro.glass PTFE gasket SS screen Silicone stopper	Boro.glass PTFE gasket SS screen Silicone stopper	Boro.glass PTFE gasket SS screen Silicone stopper


Filtration manifold Fit 250ml sterile funnel



SPECIFICATION

Item	VREZ-1EF	VREZ-3EF	VREZ-6EF
Part.No.	278105-10	278305-11	278605-12
Branch number	1 Port	3 Ports	6 Ports
Holder number	1 Set	3 Sets	6 Sets
Filter diameter	47-50mm	47-50mm	47-50mm
outlet	ø 8 mm	ø 8 mm	ø 8 mm
Parts	Material	Material	Material
Rack	SUS316	SUS316	SUS316
Handle	Anodized Alum.	Anodized Alum.	Anodized Alum.
Funnel support base	SUS316	SUS316	SUS316
Holder (Optional)	Single-pack	Single-pack	Single-pack
Funnel & Lid	PP	PP	PP
Funnel volume	250ml	250ml	250ml
Sterile	EtO	EtO	EtO

Filtration manifold Fit 250ml sterile holder



SPECIFICATION

Item	VREZ-1HF	VREZ-3HF	VREZ-6HF
Part.No.	278106-21	278306-22	278606-23
Branch number	1 Port	3 Ports	6 Ports
Holder number	1 Set	3 Sets	6 Sets
Filter diameter	47-50mm	47-50mm	47-50mm
outlet	ø 8 mm	ø 8 mm	ø 8 mm
Parts	Material	Material	Material
Rack	SUS316	SUS316	SUS316
Handle	Anodized Alum.	Anodized Alum.	Anodized Alum.
Holder base	SUS316	SUS316	SUS316
Holder (Optional)	Single-pack	Single-pack	Single-pack
Funnel & Lid	PP	PP	PP
Funnel support base	PP	PP	PP
Funnel volume	250ml	250ml	250ml
Sterile	EtO	EtO	EtO