




2ml Screw neck ND8 Vials, Standard Opening, Screw thread 8-425, 12 x 32 mm

- 8-425 screw thread caps and vials are often used in conjunction with CTC Autosampler, Shimadzu, ThermoScientific, Varian and other autosamplers.
- Clear vial, 33 Expansion borosilicate glass (Type 1, Class A)
- Amber vial, 51 Expansion borosilicate glass (Type 1, Class B)
- Vials provide graduated write-on marking spot area, measures 0.5, 1.0 and 1.5ml
- Caps are made of high quality polypropylene to exact manufacturing tolerances.
- Septa made of the highest quality materials to ensure proper function.
- Red PTFE/White Silicone Septa.
- Pre-slit Red PTFE / White Silicone Septa.

Glass Insert

Features :

- These inserts are the very best available today for HPLC, GC, LCMS & GCMS.
- Highest quality USP TYPE I, Class A, 3.3 Expansion, Borosilicate glass.

Ordering Information

2mL Screw neck ND8 Vials, Screw thread 8-425, 12x32 mm

P/N	Description	Qty
120801	2ml screw neck ND8 Vial, Clear	100/pk
120802	2ml screw neck ND8 Vial, Amber	100/pk
120803	2ml screw neck ND8 Vial, with Marking Spot, Clear	100/pk
120804	2ml screw neck ND8 Vial, with Marking Spot, Amber	100/pk
120805	2ml screw neck ND8 Vial, Clear, Type I Class A	100/pk
120806	2ml screw neck ND8 Vial, with Marking Spot, Clear, Type I Class A	100/pk

PP Screw caps – thread 8-425 for screw neck ND8 Vials, 8mm

P/N	Description	Qty
120807	8mm PP screw thread cap, Black, center hole 5.5mm with Red PTFE/White Silicone septa 8 x 1.5mm	100/pk
120808	8mm PP screw thread cap, White, centre hole, 5.5mm	100/pk
120809	8mm PP screw thread cap, Black, centre hole, 5.5mm	100/pk
120810	Septa Red PTFE/White Silicone, 8 x 1.5 mm	100/pk
120817	8mm PP screw thread cap, White, center hole 5.5mm with Red PTFE/White Silicone septa 8 x 1.5mm	100/pk

Glass Inserts - for 12 x 32mm Screw neck ND8 Vials, 8 mm

P/N	Description	Qty
120813	150ul 5x29 mm glass insert w/PP spring bottom	100/pk
120815	250ul 5x31 mm glass insert flat bottom	100/pk
120831	100ul 5x31 mm glass insert Conical bottom	100/pk

INSTRUMENT REFERENCE

Screw neck ND8 Vials, Screw Thread 8-425

The following instruments are compatible with 2.0 mL, screw neck ND8 Vials, Screw thread 8-425

Manufacturer	Instrument/Model #
A.I.M	CPS -100, 200
Alcott	738
Alltech	570
Beckman	501, 502, 507
Carlo Erba	AS800, 8000
Dani	ALS 39.80, 86.80
Dynatech	42 vial tray
Gilson	231 XL, 232 XL, 233XL, Asted XL
Gynkotech	Gina, others
Hitachi	S6651, L7200, L7250
ICI	All except LC1600
Jasco	851-AS, AS-950, LC800&900 series
LEAP Technology	CTC A200S

Manufacturer	Instrument/Model #
Perkin Elmer	Autosystem GC 420/B 4900
Pharmacia LKB	2175-010
Precision Sampling	GC111, GC311 LC241-60
Phillips	LC-XP, 4247, 4710
S.G.E.	M280D
Shimadzu	AOC-8B/9, SIL-6, AOC 14/1400, AOC-17 SIL-6B, -9A, -8A, SIL-10°, 10Ai, 10AxL
Spark Holland	SPH 125, Marathon, Promis, Triathlon, Midas

Manufacturer	Instrument/Model #
Spectra-Physics	8875, 8880
TosoHaas	TSK-6080, AS-8010 AS-8020
Unicam	4710, LC-XP, 4247
Rainin Dynamax	Dynamax AI-IA AI-200, AI-3
Varian	8100/8200, 9100/9090/9095 Marathon, 8055, 8085 9100/9090/9095/3800GC
Waters	Alliance 2690