



Screw thread ND24 Vials



Storage Vials

20mL, 30mL, 40mL and 60mL

Screw thread ND24 Vials

Borosilicate Type I Glass Vials

- These vials are ideal for use in water sampling according to EPA method.
- Clear vial, USP Type 1, Borosilicate glass.
- Amber vial, USP Type 1 Borosilicate glass.
- Caps are made from high quality polypropylene to exact manufacturing tolerances.
- Septa made of the highest quality materials to ensure proper function.
- Nature PTFE / Nature Silicone Septa.
- Blue PTFE / White Silicone Septa.
- White Silicone Septa.
- All vials are packaged in convenient trays which facilitates handling.

20mL Screw thread ND24 Vials – 27 x 58 mm ND24 Vials

P/N	Description	Qty/pk
122401	20ml screw neck ND24 vial, clear	100
122402	20ml screw neck ND24 vial, amber	100

30mL Screw thread ND24 Vials – 27 x 75 mm ND24 Vials

P/N	Description	Qty/pk
122408	30ml screw neck ND24 vial, clear	100

40mL Screw thread ND24 Vials – 27 x 95 mm ND24 Vials

P/N	Description	Qty/pk
122403	40ml screw neck ND24 vial, clear	100
122404	40ml screw neck ND24 vial, amber	100

60mL Screw thread ND24 Vials – 27 x 140 mm ND24 Vials

P/N	Description	Qty/Pk
122405	60ml screw neck ND24 vial, clear	100
122406	60ml screw neck ND24 vial, amber	100

Screw Cap and Septa Thread ND24 Vials

P/N	Description	Qty/pk
122515	24mm PP Screw Cap, white, centre hole with nature PTFE/nature Silicone septa	100
122416	24mm PP Screw Cap, white, Closed Top with nature PTFE/nature Silicone septa	100
122411	24mm PP Screw Cap, white, centre hole 15mm	100
122412	24mm PP Screw Cap, white, closed top	100
122409	Septa nature PTFE/nature Silicone, 22 x 3mm	100
122410	Septa White Silicone, 22x3mm	100