

Crimp Neck ND20 Vials



Crimp neck ND20 Vials

Borosilicate Glass USP Type I

Crimp Neck Headspace Vials are available with either a flat or rounded bottom.

- Bevelled top design.
- Available with 20ml or 10ml vial.
- Magnetic Aluminium cap suit for CTC instrument.
- Silicone/PTFE are all Available
- Rounded or flat bottom profiles are available for accurate instrument compatibility.

Crimp Neck ND20 Vials - 22.5 x 46 mm and 22.5 x 75.5 mm ND20 Vials

| Part No. | Description | Qty/pk |
|----------|--|--------|
| 122001 | 10ml crimp neck ND20 vial, Clear, flat bottom | 100 |
| 122002 | 20ml crimp neck ND20 vial, Clear, flat bottom | 100 |
| 122003 | 10ml crimp neck ND20 vial, Clear, Round bottom | 100 |
| 122004 | 20ml crimp neck ND20 vial, Clear, Round bottom | 100 |



122014

Aluminium Crimp Seals – 20 mm, for Crimp Neck ND20 Vials.

| P/N | Descriptions | Qty/pk |
|--------|--|--------|
| 122006 | Septa White PTFE/White Silicone, 20x3mm | 100 |
| 122007 | Septa Nature PTFE/Nature Silicone, 20x3mm | 100 |
| 122011 | 20mm Aluminium Crimp Cap, centre hole 10mm | 100 |
| 122012 | 20mm Blue Magnetic Aluminium Crimp Cap Cap, center hole 10mm | 100 |
| 122013 | 20mm Aluminium Crimp Cap, center hole with White PTFE/White Silicone septa | 100 |
| 122014 | 20mm Aluminium Crimp Cap, center hole with Nature PTFE/Nature Silicone septa | 100 |
| 122018 | 20mm Blue Magnetic Aluminium Crimp Cap, center hole with White PTFE/White Silicone septa | 100 |
| 122019 | 20mm Blue Magnetic Aluminium Crimp Cap, center hole with Nature PTFE/Nature Silicone septa | 100 |



20mm Hand Crimper

Crimper / Decapper – 20 mm, for Crimp Neck ND20 Vials.

| Part No. | Description | Qty/pk |
|----------|--------------------|--------|
| 122026 | 20mm Hand Crimper | 1 |
| 122027 | 20mm Hand Decapper | 1 |