

THERMOSTATIC WATER BATH



WB series thermostatic water bath adopts LCD and PID temperature control mode, with fast temperature rise and accurate temperature control. This series of thermostatic water bath pot has the function of preventing dry burning, which can avoid damage to the inner tank and heating pipe of the water bath pot caused by dry burning. Built in drainage pipe for easy drainage.

PRODUCT FEATURES

- The instrument adopts LCD display with clear and beautiful display contents, accurate and reliable, PID temperature control technology, fast temperature rise and accurate temperature control;
- The inner tank is made of high-quality stainless steel and stamped at one time. It is high-temperature resistant, well sealed and durable;
- High quality temperature sensor is adopted, with accurate temperature measurement and long service life;
- With magnetic stirring function, it can accelerate the uniform water temperature;
- Standard timing and anti dry burning function, safe and reliable.



Rapid heating



LCD Display



Pipe Weep



1-99h59
min59s



2-6 holes

Product Model	LCH-WB-2/2+	LCH-WB-4/4+	LCH-WB-6/6+
Image			
Voltage V	AC110/220		
Frequency Hz	50-60		
Power W	600	1000	1500
Display Method	LCD		
Number Of Holes	2	4	6
Temperature Control Range °C	RT~100		
Temperature Control Accuracy °C	≤±1		
Timing Range	1-99h59min59s		
Stirring Function	2+/4+/6+ have		
Pipe Weep	/	yes	yes