

# BENCHTOP LOW SPEED CENTRIFUGE



LCE series centrifuges use centrifugal force to separate liquid and solid particles or components in liquid and liquid mixtures. Separate the solid particles from the liquid in the suspension; Separate two liquids of different densities from each other in an emulsion.



## PRODUCT FEATURES

- ABS high strength shell, surface frosted, not easy to age, beautiful and light, and in line with ergonomics;
- LCD liquid crystal display, adjustable speed, time and other parameters, functional partition display, beautiful and clear ;
- DC brushless motor, low noise, maintenance free, high-speed bearing, stable operation;
- Multiple time settings, precise control, high work efficiency, multi-functional encoder adjustment, and barrier free function setting;
- With program preset function, it can save 5 different running programs ;

| Product Model              | LCH-LCF-05           | LCH-LCF-05D |
|----------------------------|----------------------|-------------|
| Image                      |                      |             |
| Voltage v                  | AC100-240            |             |
| Power W                    | 150                  |             |
| Display                    | LED                  | LCD         |
| Speed Range rpm            | 100-5000             |             |
| Speed Control Accuracy rpm | ±10                  |             |
| Timing Range               | 10s~99min59s         |             |
| Motor Type                 | DC Brushless motor   |             |
| Max RCF xg                 | 2795                 |             |
| Tube Capacity              | 8*15ml/8*12ml/8*10ml |             |