

High Speed Mini Centrifuge Model LC-Mini-10K Pro

High Speed Mini Centrifuge

Features:

- Touch button, waterproof and durable;
- Using LCD liquid crystal display, simple operation, clearer and more intuitive display;
- The use of high-strength transparent cover is convenient for real-time observation;
- Standard replaceable PCR board and angle rotor;
- Unique rotor design makes it safer to rotate and easy to replace the rotor;
- It has the functions of emergency stop when opening the lid and reversing when the lid is closed;
- The DC permanent magnet motor is durable, low-noise, stable, safe and reliable.



Specifications

Model	LC-Mini-10K Pro
Voltage	AC110-240V / 24 VDC
Frequency/Hz	50-60
Power/W	36
Display Mode	LCD
Max Speed/rpm	1000 - 10000
Max relative centrifugal force/xg	5600
Timing Range	1min - 1h 59min
Motor Type	DC permanent magnet motor
Standard Rotor Specification	8 × 2.0ml / 1.5ml / 0.5ml / 0.2ml Angle rotor 2 × 8 × 0.2ml PCR centrifugal
Net Weight/Kg	1.14
Dimensions/mm	180 × 145 × 110

