



Overhead Stirrer Model LCH-OES-30LFR

Features:

1. The LCD screen displays the set speed and actual speed, which is intuitive and clear;
2. It can realize forward and reverse stirring, stop and start slowly in the middle of forward and reverse rotation;
3. Adopt DC brushless motor, low noise, maintenance-free, stable operation;
4. Constant speed, torque can be adjusted automatically according to the change of sample viscosity;
5. Adjustable stirring rod chuck, easy to load and unload;
6. Stable start-up, effectively preventing sample overflow;
7. Automatically memorize the set working parameter.

Specification

Model	LCH-OES-30LFR
Voltage	AC100-240V, DC24V/10A
Frequency Hz	50 - 60
Power W	200
Timing range	1 - 99h59min
Display method	LCD
Speed Range rpm	50 - 1800
Max Stirring Volume (H ₂ O) L	30
Stirring Max Viscosity mPa·s	50000
Lifting Range mm	300 - 600
Motor Type	DC Brushless Motor

Ordering Information

P/N	Description
LCH-OES-30LFR	Overhead Stirrer 30L, 50-1800rpm