



SYNF10401LA

Synthware™ 47 mm filtration apparatus 1000ml with stainless steel screen support, complete joint: ST/NS 40/35



Using 47 mm membranes for filtering corrosive solutions, particulate removal from HPLC solvents, and filtration analysis with solid phase extraction disks.

Support base has removable 47 mm stainless steel support screen for use with viscous or proteinaceous solutions.

Complete assembly includes

- 300mL graduated funnel
- aluminum clamp
- Base, 47mm SS Screen Support
- 1 liter filtering flask

Inner Male Joint 40/35;

Hose Connection = 3/8" I.D. / 10mm OD

Specifications

Color	Clear
Material	Borosilicate Glass 3.3
Includes	Stainless Steel Filter
Autoclavable	Autoclavable
Graduated	Yes
Sterility	Non-sterile
Surface Treatment	Non-treated
Capacity (Metric)	1 L
For Use With (Application)	Organic Synthesis
Diameter (Metric) Disc	47 mm
Disposable	Reusable
Mouth Style	Beaded
Bottom Shape	Flat

Ordering Information

P/N	Description
F10401LA	Synthware filtration apparatus 1000ml with stainless steel screen support, complete

Replacement Part

F100001	Funnel, 47mm, 300mL Graduated
F100002	Clamp, 47mm Aluminum
F100003A	Base, 47mm SS Screen Support
F100004	Flask, 1L Filtering, 40/35
F100007	Stainless Steel Screen for Base Support

Synthware is a trademark of Kemtech America Inc.