



Product Features

1. Beautiful appearance design, and ergonomic;
2. The rubber sole is adopted to reduce the vibration during the operation of the instrument and support the smooth operation of the instrument;
3. Light weight, easy to move, save laboratory space, convenient repair and maintenance;
4. The body surface is equipped with a porous heat sink, heat dissipation effect is good, effectively improve the service life;
5. Easy to use, oil free, clean and pollution-free, simple maintenance and maintenance, can operate by itself;
6. The shell uses aluminum alloy, using electrostatic spray processing technology, durable, aging resistance, corrosion resistance.
7. The instrument uses AC asynchronous motor, long life, motor operation efficiency is outstanding and reliable;
8. The single table model is equipped with a silencer, which can reduce the noise generated during the operation of the instrument.
9. The body is equipped with fan cooling and exhaust system, which can ensure the long time operation of the instrument.

Specification

Model	LCH-VP-85 Oil free Vacuum pump
Pumping Speed (L/min)	30 l/min
Limit Pressure (mBar)	-0.085 mBar
Number of stage	1
Motor Power (Watt)	165
Dimension (L x W x H) mm	290 x 160 x 280 mm
Weight Kg	5 Kg
Noise (dB)	< 54 dB
Power	220V / 50Hz

Ordering Information

Part No. LCH-VP-85 **Oil free Vacuum pump** 30l/min, 220V/50Hz