

12-position, 24-position Water bath nitrogen blow Sample Concentrator



JTT-DC12

Specification

Product Part No.	260506	260507
Model	JTT-DC12	JTT-DC24
Sample processing capacity	12 independent control	24 independent control
Heating method	Water Bath	Water Bath
Temperature control method	PID	PID
Temperature control accuracy	± 1 °C	± 1 °C
Temperature Control range	Room Temp - 100 °C	Room Temp - 100 °C
Timer	0 - 99Hr 59min	0 - 99Hr 59min
Gas flow rate	0 – 15 l/min	0 – 15 l/min
Nitrogen consumption	330 ml/min	330 ml/min
Sample Tube size range	10 - 55mm	10 - 29mm
Sample Rack Lifting Height	0 – 150mm	0 – 150mm
Sample Rack Rotation Direction	360 °	360 °
Overall weight	10 Kg	12 Kg
Overall Dimension	420 x 545 x 675mm	420 x 545 x 675mm