

GL45 Reagent Bottle 1000ml, Borosilicate Glass 3.3



Shu Niu Original GL 45 Laboratory Bottle clear, with screw cap and pouring ring, PP, Orange , 1000 mL

- Manufactured from Shu Niu Borosilicate glass 3.3. , classified as a Glass Type I according to USP
- Standardised GL45 thread and matching cap system for particularly tight sealing and simple, clean pouring out
- Fired-on, highly durable white print
- Suitable for long-term storage
- Easy-to-read scale with graduation marks to estimate the content

Ordering Information

| Part No | Descriptions |
|---------|---|
| 260111 | Shu Nie, GL45 Reagent Bottle 1000ml, Borosilicate 3.3 |