



P/N 250201 GL45 Safety Waste cap with activated carbon filter module and HDPE waste tank 5 liters



P/N 250202 Activated carbon filter module With 6 month time strip



Application Scenarios

HPLC waste collection/AAS waste collection/LC-MS waste collection

Product Features


1. G45 universal threaded connector to match conventional laboratory waste drums.
2. High-efficiency filter, effectively absorbing toxic and harmful gases, balancing the pressure of the waste liquid tank.
3. Multi-functional waste liquid piping connector can match different capillary tubes and silicone hoses.
4. Waste liquid is sealed and recycled, making your lab safer.

Purchase Instructions:

1. The standard configuration of the 1 x PP & PTFE, G45 thread, 1×M8 barbed connector port, 2 x 1/4-28 PP thread adapter for tubing OD 1/8", 1 x Exhaust filter, activated carbon with 6 month time strip
2. The G45 Safety Waste Cap can be used with Waste bottle GL45 thread.

NOTICE:

1. When you first start using it, just press your finger on the time bar to start it;
2. The time bar will gradually turn red after startup. When it turns completely red, it is time to replace the Exhaust Activated Carbon module, as shown in the figure below;
3. We recommend replacing the Activated Carbon filter module after 6 months of use.



Ordering Information

Part No.	Descriptions	Qty
250201	G45 Safety waste cap completed with : 1 x PP & PTFE GL45 thread, 1 x M8 barbed connector port, 2 x 1/4-28 PP thread adapter for tubing OD 1/8", 1 x Exhaust activated carbon filter module with 6 month time strip 1 x HDPE Wasted tank 5 liter	1

Accessories

250202	Exhaust Activated Carbon filter module with 6 month time strip	1
0522	HDPE Waste tank 5 liter, H x W x D 24.0 x 14.5 x 21.4 mm	1