

GL45 HPLC Mobile Phase Kit, 2 port, With 1000ml Reagent Bottle, completed



Features & Benefits

- Convenient kits for easy mobile phase management.
- Protect your mobile phase from evaporation and contamination, and your lab from harmful vapors.
- **P/N 260100** GL45 HPLC Mobile Phase Kit, 2 port complete with :
 - 1 x GL45 Borosilicate Glass 3.3 Reagent bottle 1000ml
 - 1 x GL45 Mobile Phase Cap 1/4-28 thread 2 port
 - 2 x PEEK Connector 1/4-28 thread for PTFE Tubing 1/8" OD
 - 1 x 3 meters PTFE Tubing OD 1/8"
 - 1 x 10 µm SUS 316 Filter
- Use one 1/4-28 threaded port for exhaust filter or air inlet valve; other ports can be used as needed or closed with 1/4-28 plug (sold separately).

Ordering Information

Part No	Descriptions
260100	GL45 Safety Cap for HPLC Mobile Phase Kit, 2 port, With 1000ml Reagent Bottle, completed

Spare Part

Part No	Descriptions
260111	Reagent Bottle 1000ml, Borosilicate 3.3
260112	GL45 Cap, 1/4-28 thread 1 port
260113	GL45 Cap, 1/4-28 thread 2 port
260114	GL45 Cap, 1/4-28 thread 3 port
260115	GL45 Cap, 1/4-28 thread 4 port
260116	PEEK Connector 1/4-28 thread for Tubing 1/8"OD.
260117	PEEK Seal Plug Connector 1/4-28 thread
260118	PP Seal Plug Connector 1/4-28 thread
260119	One way valve with inbuilt PTFE membrane
260120	SUS 316L Filter head for Tubing 1/8" OD.
260121	Glass Filter head for Tubing 1/8" OD.
260122	PP Ion Chromatograph Filter head for Tubing 1/8" OD.
260123	Time Strip 6 month
260124	3 meter PTFE Tubing OD 1/8"