

## GL45 HPLC Mobile Phase Cap, 1/4-28 thread 1, 2, 3 port, with PEEK Connector 1/4-28 thread for Tubing 1/8" OD.



260125



260126



260127

### Features & Benefits

- Protect your mobile phase from evaporation and contamination, and your lab from harmful vapors.
- Includes PEEK Connector 1/4-28 thread nuts with ferrules.
- Use one 1/4-28 threaded port for exhaust filter or air inlet valve; other ports can be used as needed or closed with 1/4-28 plug (sold separately).

### Ordering Information

Part No	Descriptions
260125	GL45 Cap, 1/4-28 thread 1 port with PEEK Connector for Tubing 1/8" OD
260126	GL45 Cap, 1/4-28 thread 2 port with PEEK Connector for Tubing 1/8" OD
260127	GL45 Cap, 1/4-28 thread 3 port with PEEK Connector for Tubing 1/8" OD
260116	PEEK Connector 1/4-28 thread for Tubing 1/8"OD.
260117	PEEK Seal Plug Connector 1/4-28 thread



260116



260117