

## GL45 Safety Cap for HPLC Mobile Phase



### Features & Benefits

- Protect your mobile phase from evaporation and contamination, and your lab from harmful vapors.
- Includes GL45 cap 1/4"-28 thread 2 port, one way valve with inbuilt PTFE membrane and PEEK Connector 1/4"-28 thread nuts with ferrules.

### Ordering Information

Part No	Descriptions
260103	GL45 Safety Cap for HPLC Mobile Phase Kit, 2 port, With one way valve and Time strip

### Spare Parts

260112	GL45 Cap, 1/4"-28 thread 1 port
260113	GL45 Cap, 1/4"-28 thread 2 port
260114	GL45 Cap, 1/4"-28 thread 3 port
260115	GL45 Cap, 1/4"-28 thread 4 port
260116	PEEK Connector 1/4"-28 thread for Tubing 1/8"OD.
260117	PEEK Seal Plug Connector 1/4"-28 thread
260118	PP Seal Plug Connector 1/4"-28 thread
260119	One way valve with inbuilt PTFE membrane
260120	Time Strip 6 month



260119



260120



260113