

Ultrasonic Cleaner Model KZ-F6



Specification

Model	KZ-F6
Capacity	6.5 Liter
Frequency	40 KHz
Transducer	3 pcs
Timer	0 - 30 minutes
Degasser	Yes
Temperature Range	+5 °C room temp to 80 °C
Heater (Watts)	300 watts
Ultrasonic (Watts)	180 watts
Inner Tank Dimension L x W x H mm	300 x 150 x 150 mm
Overall Dimension L x W x H mm	395 x 225 x 325 mm
Tank Material	SS 304
Basket Material	SS 304
Weight Kg	4.6 Kg
Power	220 V / 50Hz