



CeiExpert Boro 3.3 Filtration Set with Flask 1000ml

The vacuum filter is a laboratory equipment, in the technical chromatographic, for the removal of column life and impurity detection system accuracy, plays an important role. In the filtration process, the liquid degassing effect is obvious. The equipment for different purposes of membrane filtration has been applied in weight analysis, chromatography, colloidal separation and other tests.

Main parameters

Material	: Borosilicate glass 3.3
Cup capacity	: 300 ml
Flask Bottle receiving capacity	: 1000 ml
Connector	: 40 mm
Outer diameter of sand core	: 60 mm
Inner diameter of sand core	: 42 mm
Use with Membrane filter diameter	: 47 mm / 50 mm
Tubing	: Silicone Tubing 7mm ID x 14mm OD 1 meter

Ordering Information

Part No	Descriptions
260201	CeiExpert Boro 3.3 Filtration set with Flask 1000mL
260202	CeiExpert Boro 3.3 Filtration set with Flask 2000mL
260205	CeiExpert , 300mL Funnel
260206	CeiExpert , Flask bottle 1000ml
260207	CeiExpert , Sand Core Filter Head, 40/38
260208	CeiExpert , Aluminum Clamp
260209	Silicone Tubing 7mm. ID x 14 mm OD 1 meter