



Glass insert

2.0 mL, 9 mm Wide Opening, Short Cap, Screw-Thread Vials 12 x 32 mm

Features :

- Wider opening prevents syringe problems.
- 9mm screw thread caps, septa, and vials are designed especially for use with Agilent Technologies, AB Sciex, Brukers, CTC PAL, Dionex, Hitachi, Knauer, PerkinElmer, Shimadzu, ThermoScientifics, Waters, Alliance's and other rotating or robotic arm samplers.
- Clear vial, 33 Expansion borosilicate glass (USP Type 1, Class A)
- Amber vial, 51 Expansion borosilicate glass (Type 1, Class B)
- Vials provide graduated write-on marking spot area, measures 0.5, 1.0 and 1.5ml
- Caps are made of high quality polypropylene to exact manufacturing tolerances.
- Septa made of the highest quality materials to ensure proper function.
- Red PTFE/White Silicone Septa.
- Pre-slit Red PTFE / White Silicone Septa.

Glass Insert

Features :

- These inserts are the very best available today for HPLC, GC, LCMS & GCMS.
- Highest quality USP TYPE 1, Class A, 33 Expansion, Borosilicate glass.

Ordering Information :

9mm Vial ND9, Wide opening with Short Thread

120901	2ml screw neck ND9 vial, Clear	100/pk
120902	2ml screw neck ND9 vial, Amber	100/pk
120903	2ml screw neck ND9 vial with Marking spot, Clear	100/pk
120904	2ml screw neck ND9 vial with Marking spot, Amber	100/pk
120905	2ml screw neck ND9 vial, Clear, with integrated 300 ul Micro glass insert, based bonded	100/pk

9mm Glass Insert

120913	250ul 5.7 x 29 mm glass insert w/PP bottom spring	100/pk
120915	300ul 6 x 31 mm glass insert flat bottom	100/pk
120925	150ul, 5.7 x 29 mm glass insert conical bottom	100/pk

9mm Screw cap and Septa

120907	9mm PP screw thread cap, Blue center hole with Red PTFE/White Silicone septa	100/pk
120908	9mm PP screw thread cap, Blue center hole with White PTFE/Red Silicone septa	100/pk
120909	9mm PP screw thread cap centre hole, 6mm	100/pk
120910	Septa Red PTFE/White Silicone, 9 x 1mm	100/pk
120911	Pre-Slit Septa Red PTFE/White Silicone, 9 x 1mm	100/pk
120912	Septa White PTFE/Red Silicone, 9 x 1mm	100/pk
120914	Pre-Slit Septa White PTFE/Red Silicone pre slit, 9 x 1mm	100/pk
120916	9mm PP screw thread cap, Close-top	100/pk
120917	9mm PP screw thread cap, Blue center hole with Red PTFE/White Silicone/Red PTFE septa	100/pk
120919	9mm PP screw thread cap, Blue center hole with Pre-Slit White PTFE/Red Silicone septa	100/pk
120920	9mm PP screw thread cap, Blue center hole with Pre-Slit Red PTFE/White Silicone septa	100/pk
120922	Septa Red PTFE/White Silicone/Red PTFE, 9 x 1mm	100/pk

Rack for vial 2ml

120109	Vial rack for 2 ml vials, Blue, 50 holes	1/pk
--------	--	------

INSTRUMENT REFERENCE

For 2.0 mL, 9mm Screw neck ND9 Vials

The following instruments are compatible with 2.0 mL, 9mm screw neck ND9 Vials

Manufacturer	Instrument/Model #
Carlo Erb/Fison	AS-800
Dani	ALS 39.80, 86.80
Dynatech	42 vial tray, LC-200
Gilson	231 XL, 232 XL, 233XL, Aspec XL
Agilent Technologies	
Hitachi	AS-2000, AS-6000
LCD	Marathon, Promis
L.E.A.P. Technologies	CTC A200S CTC A105S

Manufacturer	Instrument/Model #
Perkin Elmer	AI-1, Integral 4000 ISS100, ISS200, LC600 42 vial tray
Shimadzu	SIL-10A, 10Ai, 10AxL
Spark Holland	SPH 125, Marathon, Promis
Spectra-Physics	8875, 8880
Thermo Scientific	TRACE GC 2000 AS2000
Varian	8100/8200, Marathon
Waters	Alliance 2690