



Specification

Model	E1.3L	E2L	E3.2L	E4.5L	E6.5L	E10L	E15L	E22L	E30L
Capacity liter	1.3	2	3.2	4.5	6.5	10	15	22	30
Ultrasonic Power Watts	60	60	120	180	180	240	360	480	600
Heater Power Watts	100	100	150	300	300	300	450	600	600
No. of Transducer	1	1	2	3	3	4	6	8	10
Frequency Hz	40								
Temperature Range °C	+5 °C from room temp – 80 °C								
Timer min	0 – 30 min								
Degasser	Yes								
Inner Dimeter mm	150x135x60	150x135x100	240x135x100	300x150x100	300x150x150	300x325x150	325x300x155	500x300x150	500x300x200
Weight Kg	1.9	1.9	2.95	4.2	4.6	6	8.4	10.35	12
Outer Dimeter mm	175x160x165	175x160x195	265x165x225	325x175x235	325x175x275	330x270x275	355x325x275	530x330x275	530x330x330

Ordering Information

P/N	Description	ราคา บาท***
260901	Model E1.3L Ultrasonic cleaner 1.3 liters	4,500
260902	Model E2L Ultrasonic cleaner 2 liters	4,500
260903	Model E3.2L Ultrasonic cleaner 3.2 liters	6,300
260904	Model E4.5L Ultrasonic cleaner 4.5 liters	9,000
260905	Model E6.5L Ultrasonic cleaner 6.5 liters	11,000
260906	Model E10L Ultrasonic cleaner 10 liters	12,000
260907	Model E15L Ultrasonic cleaner 15 liters	14,500
260908	Model E22L Ultrasonic cleaner 22 liters	18,000
260909	Model E30L Ultrasonic cleaner 30 liters	23,000

*** ราคาไม่รวมภาษีมูลค่าเพิ่ม 7%