



Specification

Model	HH-1	HH-2	HH-4	HH-6	HH-8
Power	220V/50Hz				
Display	LED				
Temperature Range	+5 °C from room temp – 100 °C				
Temperature Accuracy	± 1 °C				
Timer	999 min				
Position	1	2	4	6	8
Inner Dimeter mm	160x160x120	300x160x120	300x300x120	470x300x120	630x320x120
Weight Kg	3.2	4.6	6.8	8.6	7.4
Outer Dimeter mm	200x200x200	350x200x200	350x350x200	510x350x200	740x360x180

Ordering Information

P/N	Description	ราคา บาท***
260801	Model HH-1 Water Bath 1 hole	3,500
260802	Model HH-2 Water Bath 2 hole	4,500
260803	Model HH-4 Water Bath 4 hole	6,000
260804	Model HH-6 Water Bath 6 hole	7,500
260805	Model HH-8 Water Bath 8 hole	12,000
260806	Model HH-1 Water Bath 1 hole with anti dry burning	4,800
260807	Model HH-2 Water Bath 2 hole with anti dry burning	5,800
260808	Model HH-4 Water Bath 4 hole with anti dry burning	7,500
260809	Model HH-6 Water Bath 6 hole with anti dry burning	11,000

*** ราคาไม่รวมภาษีมูลค่าเพิ่ม 7%